



Living up to Life

Leica
MICROSYSTEMS

Material
Research





Life Science
Research

Application Note


Targeting of peroxisomal matrix proteins in the diatom
Phaeodactylum tricornerutum

related instrument Leica EM PACT2 / EM AFS2

Medical
Research



Industrial
Manufacturing



Natural
Resources



Targeting of peroxisomal matrix proteins in the diatom *Phaeodactylum tricornerutum*

COURTESY


Dr. Kathrin Bolte and Dr. Andreas Klingl
Zellbiologie & Elektronenmikroskopie
AG Prof. Uwe-G. Maier
Philipps Universität Marburg, Marburg
The website of the EM facility is:
<http://www.synmikro.com/de/startseite/core-facilities/elektronenmikroskopie.html>

SAMPLE PREPARATION:

P. tricornerutum cells expressing different types of GFP fusion proteins were harvested via centrifugation at 1,500xg and cryoimmobilized by high-pressure freezing on gold carriers (Leica EM PACT2, Leica Microsystems GmbH, Vienna, Austria). Subsequently, the cells were freeze-substituted with acetone in combination with 2% OsO₄, 0.25% uranyl acetate and 5% H₂O. Freeze substitution was carried out in the automated Leica EM AFS2 unit (Leica Microsystems GmbH, Vienna, Austria) at -90°C for 4 h, -60°C for 8 h, -30°C for 8 h and then held at 0°C for at least 3 h. The heating time between each step was 1 h. After washing the samples in ice-cold acetone, they were gradually infiltrated in Epon at room temperature, followed by polymerization at 60°C for three days. Ultrathin sections of embedded samples were collected on uncoated nickel grids (400 square mesh). Transmission electron micrographs were either taken on a JEOL 2100 TEM operated at 80 kV in combination with a fast-scan 2K x 2K CCD camera F214 (TVIPS, Gauting, Germany) or on a Zeiss CEM 902 operated at 80 kV equipped with a wide-angle Dual Speed 2K CCD camera (TRS, Moorenweis, Germany).

For more details see:

A single peroxisomal targeting signal mediates matrix protein import in diatoms. Gonzalez NH, Felsner G, Schramm FD, Klingl A, Maier UG, Bolte K, (2011) PLoS ONE 6(9): e25316. doi:10.1371/journal.pone

**LEGEND FOR IMAGES:**

(A) Ultrathin section of *P. tricornutum* expressing Pex10-GFP in Epon without antibody labeling. The boxed area is shown in (B) at higher magnification and illustrates two peroxisomes in proximity to the nucleus, the golgi and the plastid.

CW, cell wall; G, golgi apparatus; Mt, mitochondrium; Ne, nuclear envelope; Nu, nucleus; P, peroxisome; PI, plastid; V, vacuole. Scale bars represent 2 μm (A), 500 nm (B).

RELATED PRODUCTS



Leica EM PACT2



Leica EM AFS2

The statement by Ernst Leitz in 1907, *“With the User, For the User,”* describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

LIFE SCIENCE DIVISION - NANO TECHNOLOGY LNT

The Leica Microsystems Nano Technology Division’s focus is to provide the most comprehensive product portfolio for the preparation of biological, medical and industrial samples for investigation in the Electron and Light Microscope. Excellent Sample Preparation is a prerequisite for perfect microscopy. **Your image starts here!**

Leica Microsystems – an international company with a strong network of worldwide customer services:

Active worldwide	Tel.	Fax
Australia · North Ryde	+61 2 8870 3500	2 9878 1055
Austria · Vienna	+43 1 486 80 50 0	1 486 80 50 30
Belgium · Diegem	+32 2 790 98 50	2 790 98 68
Canada · Concord/Ontario	+1 800 248 0123	847 405 0164
Denmark · Ballerup	+45 4454 0101	4454 0111
France · Nanterre Cedex	+33 811 000 664	1 56 05 23 23
Germany · Wetzlar	+49 64 41 29 40 00	64 41 29 41 55
Italy · Milan	+39 02 574 861	02 574 03392
Japan · Tokyo	+81 3 5421 2800	3 5421 2896
Korea · Seoul	+82 2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31 70 4132 100	70 4132 109
People’s Rep. of China · Hong Kong	+852 2564 6699	2564 4163
· Shanghai	+86 21 6387 6606	21 6387 6698
Portugal · Lisbon	+351 21 388 9112	21 385 4668
Singapore	+65 6550 5999	6773 0628
Spain · Barcelona	+34 93 494 95 30	93 494 95 32
Sweden · Kista	+46 8 625 45 45	8 625 45 10
Switzerland · Heerbrugg	+41 71 726 34 34	71 726 34 44
United Kingdom · Milton Keynes	+44 800 298 2344	1908 246312
USA · Buffalo Grove/Illinois	+1 800 248 0123	847 405 0164