PEARL



The Compact Solution to **Quality Tissue Processing**

feica

HistoCore PEARL Tissue Processor

HistoCore PEARL Tissue Processor



Ultimate Protection for User Safety

- · Fully enclosed system
- Triple Fume Protection: extraction at the retort, an external exhaust adapter and an activated carbon filter provide added safety for lab staff

Ultra-Compact Solution

- 200 cassette capacity
- Small footprint to maximize laboratory space
- 3.5L reagent bottles and three paraffin baths help maximize run efficiency

Simple Operation

- Easy-to-remember operating procedures on xylene resistant, multi-lingual touch-screen interface
- Three pre-validated protocols and Reagent Management System (RMS) provides added flexibility
- USB connection to download log files help maintain accurate records

Easy-to-Clean Paraffin and Reagent Exchanges

- Drawer style design for paraffin baths allows for easy filling and draining
- Trays under the paraffin baths and reagent bottles effectively collect spilled liquids/chemicals

















LeicaBiosystems.com





Contact your Leica Biosystems representative for a full range of validated consumables for the HistoCore PEARL Tissue Processor.

SYSTEM SPECIFICATIONS

GENERAL		Capacity	Max. 3.5 L per paraffin bath
Nominal supply voltages	100 to 120 V AC or 220 to	Melting time	50 °C to 64 °C, 12 hours
	240 V AC	Temperature setting range	50 °C to 70 °C
Nominal frequency Nominal power	50 / 60 Hz 1650 VA	BOTTLES	
Operating temperature range	+10 °C to +35 °C	Number of reagent bottles	11 (in reagent cabinet)
RETORT	Max. 200 cassettes	Condensate bottle Cleaning solution bottles	1 2
Capacity Temperature (paraffin)	50 °C to 70 °C	Bottle volume	Max. 3.5 L
Temperature (processing reagents)	Ambient temperature or 35 °C to 60 °C	DIMENSIONS AND WEIGHTS	
		Dimensions, unpacked	
Temperature (cleaning reagents)	Alcohol 62 °C, Xylene 67 °C	(W x D x H)	600 x 710 x 1320 mm (23.62 x 27.95 x 51.96 in)
PARAFFIN BATHS		Weight, unpacked	Max. 140 kg (308.65 lbs)
Number of paraffin baths	3	Weight, packed	Max. 175 kg (385.81 lbs)

Leica Biosystems—an international company with a strong network of worldwide customer service.

North America Sales and Customer Support		
North America	800 248 0123	
Asia/Pacific Sales and Customer Support		
Australia	1800 625 286	
China	+85 2 2564 6699	
Japan	+81 03 6758 5690	
South Korea	+82 2 3416 4500	
New Zealand	0800 400 589	
Singapore	+65 6550 5999	
Latin America Sales and Customer Support		
Brazil	+55 11 2764 2411	
Europe Sales and Customer Support		

For detailed contact information about European sales offices or distributors, please visit our website: LeicaBiosystems.com

Contact your Leica Biosystems representative today to learn more about the HistoCore PEARL Tissue Processor.

LEICA BIOSYSTEMS

 $\label{leader} Leica\ Biosystems\ is\ a\ global\ leader\ in\ workflow\ solutions\ and\ automation, integrating\ each\ step\ in\ the\ workflow\ from\ biopsy\ to\ diagnosis.$



Copyright © 2016 by Leica Biosystems, Nussloch, Germany. All rights reserved. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH.