

ET-101/1 GROSSING STATION / 1 USERS



The grossing stations ET-101/1 Fully stainless steel construction. Absorption system maintaining a clean atmosphere in the workplace. Dimensions: length 1300 mm, depth 850 mm, height 2100 mm. Made-up of: front panel, working surface, suction system *and* bottom panel.

Features:

- Fully stainless steel construction.
- Absorption system maintaining a clean atmosphere in the workplace.
- Air aspiration Speed 0,3 – 0,4 m³/sec
- Divided into 2 parts.
- 1 working areas 600 x 500 mm | 1 sink 400 x 400 x 250 mm. .
- Formaldehyde automatic dosing system, by pedal .
- The double suction system is situated in the frontal part and under the working area, which stimulates the suction of the contaminated air through its perforated lids (formaldehyde is heavier than air) through nozzles into regenerators **aluminum oxide filters impregnated with potassium permanganate.**
- Magnetic board for cutting Tools
- Working area with removable perforated cover. Below this working area are the suction system and the general drainage system. | Carving board 300 x 400 mm

Do you need more information ? Contact us

DMI MEDICAL INC 4611 S University DR Suite #435 FL3328 | 954-538-3808