

ET-101/2 GROSSING STATION / 2 USERS



The grossing stations ET-101/2
Fully stainless steel construction. Absorption system maintaining a clean atmosphere in the workplace. Dimensions: length 2150 mm, depth 850 mm, height 2100 mm. Made-up of: front panel, working surface, suction system *and* bottom panel. General lighting system. 2 adjustable arms for illumination

Features:

- Fully stainless steel construction.
- Absorption system maintaining a clean atmosphere in the workplace.
- Air aspiration Speed 0,3 – 0,4 m³/sec
- Divided into 3 parts.
- 2 working areas 600 x 500 mm | 1 sink 400 x 400 x 250 mm. .
- Formaldehyde automatic dosing system, by pedal .
- Residual formaldehyde container with level detector alarm (capacity up to 10 liters)
- Magnetic board for cutting Tools
- Working area with removable perforated cover. Below this working area are the suction system and the general drainage system. | Carving board 300 x 400 mm

Do you need more information ? Contact us

DMI MEDICAL INC 4611 S University DR Suite #435 FL3328 | 954-538-3808