

ET-107-2 NEW GROSSING STATION, 2 USERS



The grossing stations

ET-100 SERIES form a family of products characterized by their flexibility of adaptation to space and the versatility of configuration. Special attention has been paid to the conception of the work place, in order to make the tasks of the laboratory more comfortable and safe. This allows the user to execute their work more efficiently.

Features:

- Manufactured entirely in stainless steel AISI-304, with sheet thickness between 1.5 and 2 mm.
- Front safety glass screen, with guillotine electric system
- Aspiration capacity: 2450 m³ / hour (depending on model).
- Formalin dispenser tap.
- Formalin collection sink (includes pre-filter and lid).
- Integrated waste collection bucket.
- Carving area w/removable perforated trays, which guarantee the descending aspiration in the work area.
- Auxiliary electrical outlets.
- **POSSIBILITY OF ADAPTATION TO SPECIFIC NEEDS**

Do you need more information ? Contact us

DMI MEDICAL INC 4611 S University DR Suite #435 FL3328 | 954-538-3808